


Temporary Collagen Inserts



ASSESS THE EFFECTIVENESS OF PERMANENT PUNCTAL OCCLUSION

EagleVision dissolvable punctal/canalicular Collagen Inserts can be used as a cost effective way to evaluate the benefits of permanent punctal occlusion. Collagen may be inserted to aid in diagnosis of contact lens intolerance, enhance ocular medications, and increase patient comfort after surgery.



- Easy to insert
- Expands in the punctum after insertion to block tear drainage
- Effective for 3-5 days - dissolves completely in 7-10 days
- Available in 1.6mm "shorts" and 2.0mm lengths; 3 diameters
- Helps retain ocular medications after surgery and foreign body removal
- Instructional video available

Supplied	Description	REF
6 inserts per sterile foam pack, 10 packs per box. Total of 60 inserts.	Temporary Collagen Inserts	0025 - 0.2mm diameter x 2.0mm length 0026 - 0.3mm diameter x 2.0mm length 0027 - 0.4mm diameter x 2.0mm length 0051 - Variety Pack
Variety Pack - Sterile foam packs of (4)each 0.2mm, (4)each 0.3mm, and (2)each 0.4mm diameter. Total of 60 inserts.	Temporary Collagen Inserts Variety Pack	
	Temporary Collagen "Shorts" Inserts	0052 - 0.2mm diameter x 1.6mm length 0053 - 0.3mm diameter x 1.6mm length 0054 - 0.4mm diameter x 1.6mm length 0057 - Collagen "Shorts" Variety Pack
	Temporary Collagen "Shorts" Inserts Variety Pack	



8500 Wolf Lake Dr., Suite 110 • P.O. Box 34877 • Memphis, TN 38133 USA
 1-800-222-7584 • (901) 380-7000 • Fax (901) 380-7001
www.eaglevis.com • info@eaglevis.com