



# Gellansert®

## Temporary Intra Canalicular Inserts



### ASSESS THE EFFECTIVENESS OF PERMANENT PUNCTAL OCCLUSION

EagleVision dissolvable intra canalicular Gellan Inserts can be used as a cost effective way to evaluate the benefits of permanent punctal occlusion. Gellan may be inserted to aid in diagnosis of contact lens intolerance, enhance ocular medications, and increase patient comfort after surgery. Expands to fit up to 1.2mm diameter. Easily removed with saline irrigation.



CE Marked

- Easy to insert
- Expands in the canaliculus after insertion to block tear drainage
- Intended to temporarily block tear drainage by occlusion of the canaliculus
- One size fits all
- Helps retain ocular medications after surgery and foreign body removal
- All natural – bovine free
- Made of high molecular weight polysaccharide gum
- Instructional video available

Supplied	Description	REF
5 inserts per sterile foam pack, 20 packs per box. Total of 100 inserts.	Temporary Gellan Inserts 0.3mm dia. x 1.5mm length	0040



8500 Wolf Lake Dr., Suite 110 • P.O. Box 34877 • Memphis, TN 38133 USA  
1-800-222-7584 • (901) 380-7000 • Fax (901) 380-7001  
www.eaglevis.com • info@eaglevis.com